

SEPTANIOS MD

Pre-disinfectant detergent
of instruments



- **Neutral pH: compatible with all types of materials**
- **Compatible with corrosion-resistant metals, glass, and plastics (polycarbonates, ...)**
- **Broad antimicrobial spectrum in 15 minutes**
- **Concentration for use: 0.5%**
- **Easily biodegradable oleochemical detergent**


INDICATIONS

- Cleaning and pre-disinfection of surgical and medical instruments, medical devices, heat-sensitive instruments, and endoscopic equipment.
- Collection of soiled instruments.

CHARACTERISTICS

- Clear, blue solution
- Non corrosive to instruments
- Suitable for use in ultrasonic tanks
- pH of diluted product at 0.5% : about 7


**PRODUCT
TO DILUTE**
25 ml > 5L



SEPTANIOS MD


Pre-disinfectant detergent of instruments

INSTRUCTIONS FOR USE



0.5 % dilution:


Pour a dose of 25 ml into 5 litres of cold water.
Renew the bath at least once a day.



Completely immerse the medical device.

Advised soaking time: 15 min.

Brush if necessary.
For endoscopic equipment: brush channels.



Thoroughly rinse the medical device with tap water (good microbiological grade).
For endoscopic equipment: outside and inside.



Wipe using a clean drape.
Go to the next step (see the in-house protocol).

COMPOSITION

Didecyltrimethylammonium chloride, Chlorhexidine digluconate, polyhexamethylene biguanide, easily biodegradable oleochemical non-ionic surfactant, fragrance and colorant.

PRECAUTIONS FOR USE

Dangerous. Observe the precautions for use (according to European Directive 99/45/EC and its adaptations).
Storage: between +5°C and +35°C.
Class IIb medical device (Directive 93/42/EEC as amended)

PACKAGING

- 1 12 dosing bottles of 1L.....Ref. 1915.095
- 2 4 x 5L cans + 25ml pump.....Ref. 1915.036
- 200 x 25ml doses.....Ref. 1915.097

MICROBIOLOGICAL PROPERTIES

Recommendations of use: 0,5 % - 15 min.
Adjust these conditions according to the required antimicrobial activity.

Active against	Standards	Contact time	
		0,5%	2%
Bacteria	EN 1040, EN 13727	5 min.	
	EN 14561	15 min.	
	BMR (EN 13727)	5 min.	
Mycobacteria	Mycobacterium terrae (EN 14563)	240 min.	30 min.
Yeasts	EN 1275, EN 13624, EN 1650	5 min.	
	EN 14562	30 min.	10 min.
Moulds	EN 1275	15 min.	
	EN 14562	60 min.	10 min.
Viruses	HIV-1, PRV (surrogate of HBV), BVDV (surrogate of HCV), Herpesvirus,	5 min.	
	Influenza A virus H ₁ N ₁	15 min.	
	Vaccinia virus		15 min.

