

SURFANIOS Premium

Detergent disinfectant
for floors and surfaces



- **Twin action: combination of detergent and disinfectant action in a simultaneous operation**
- **Activity against bacteria and moulds isolated from hospital environment**
- **Aldehyde free product**
- **Optimisation of ecotoxicological information**

INDICATIONS

Cleaning and disinfection of floors, walls, equipment including medical equipment.

CHARACTERISTICS

- Green coloured solution
- pH of pure product : 11.2
- pH of diluted product (0.25%) approx 8
- Density: 1.007
- Refraction index: 1.3587

PRODUCT
TO BE DILUTED
20 ml > 8 L



SURFANIOS Premium

Detergent disinfectant for floors and surfaces


INSTRUCTION FOR USE



1 Fill a washing bucket with 8 litres of water. **0.25 %** dilution: Pour a 20 ml dose of SURFANIOS PREMIUM in the washing bucket.



2 Proceed to clean respecting the draw: from the back to the exit (do not rinse surfaces).



3 Change the mop head as often as necessary. Place the dirty mop head in the dedicated container.



4 Replace the mop head after cleaning each area and change the washing solution.

COMPOSITION

N-[3-aminopropyl]-N-dodecylpropane-1,3-diamine – CAS N° 2372-82-9 (51 mg/g), didecyldimethylammonium chloride – CAS N° 7173-51-5 (25 mg/g), excipients.

PRECAUTIONS FOR USE

Dangerous. Follow the instructions for use (in accordance with Directive 99/45/EC and its adapted versions). Use the biocides with precaution. Before use, read the label and the product information. Store between +5°C and +35°C. Biocidal product for disinfecting surfaces (Main group 1 – PT 2 & 4 - Directive 98/8/EC). Reserved for professional use.

PACKAGING

SURFANIOS Premium

- 1 500 dosis, 20 ml each.....Ref. 1917.129
- 2 12 dosing bottles, 1L each.....Ref. 1917.092
- 3 4 cans of 5 litres with 1 dosing pump of 20 ml.....Ref. 1917.036

MICROBIOLOGICAL PROPERTIES

Active against	Standards	Contact time		
		0,25%	0,5%	1%
Bacteria	EN 1040, EN 1276, EN 13727	5 min.		
	BMR (EN 1276)	5 min.		
	EN 13697	30 min.	15 min.	
	BMR (EN 13697)	15 min.		
Mycobacteria	Mycobacterium terrae (EN 14348)	120 min.		
	Mycobacterium terrae (EN 14563)	120 min.	60 min.	
Yeasts	EN 1275, EN 1650	15 min.		
	EN 13697	30 min.	15 min.	
Moulds	EN 1275	15 min.		
	EN 13697	30 min.		15 min.
	EN 1650	120 min.		
Viruses	According to EN 14476 HIV-1, PRV (surrogate of HBV), BVDV (surrogate of HCV) and Herpesvirus	5 min.		
	Vacciniavirus	60 min.	15 min.	

