SURFANIOS SURFANIOS CITRON



- Twin action: combination of detergent and disinfectant action in a simultaneous operation
- Activity against bacteria and moulds isolated from hospital environment
- Active against Legionella pneumophila
- Aldehyde free formulation



INDICATIONS

Cleaning and disinfection of floors, walls and medical equipment.

CHARACTERISTICS

- Clear blue-green coloured solution with perfume.
- Usable in cold or hot water (up to +60°C).
- Wide compatibility with materials and surface coatings: pH of diluted product near to neutrality.
- pH after dilution (0.25%) with tap water: approximately 8.5
- Non corrosive (without oxydising agents).







SURFANIOS CITRON

Detergent disinfectant for floors and surfaces

INSTRUCTION FOR USE



Fill a washing bucket and a rinsing bucket with 8 litres of water. **0.25 % dilution:** Pour a 20 ml dose of SURFANIOS in the washing bucket.



Proceed with cleaning respecting the draw: from the back to the exit. Rinsing is not necessary unless the surface comes into contact with skin or mucosa.



Rinse and wring the mop before re-soaking in the washing bucket. (the cleaning trolley has to stay in the corridor, outside of the room).

COMPOSITION

N-(3-aminopropyl)-N-dodécylpropane-1,3-diamine – CAS N° 2372-82-9 (51 mg/g), didecyldimethylammonium chloride- CAS N° 7173-51-5 (25 mg/g), excipients.



PRECAUTIONS FOR USE

Dangerous. Follow the instructions for use (Drawn up according to the European rules in force regarding the classification and labelling of chemical products) Use the biocides with precaution. Before use, read the label and the product information. Store between +5°C and +35°C. Biocidal product intended for disinfection of surfaces (Group 1 - PT 2) - for professional use only.



PACKAGING

SURFANIOS

1 500 dosis, 20 ml each	.Ref.	350.129
2 12 dosing bottles, 1L each	Ref.	350.092
3 4 cans of 5 litres with 1 dosing pump of 20 ml	.Ref.	350.036

SURFANIOS CITRON

500 dosis, 20 ml each	Ref. 347.129
12 dosing bottles, 1L each	Ref. 347.092
4 cans of 5 litres with 1 dosing pump	of 20 mlRef. 347.036

MICROBIOLOGICAL PROPERTIES

Active against	Standards	Contact time
Bacteria	EN 1040, EN 13727, EN 1276 T 72-300 (MdRB)	5 minutes
	T 72-300 (L. pneumophila)	15 minutes
Mycobacteria	M. terrae (EN 14348, EN 14563)	30 minutes / 0.5%
	M. terrae (EN 14348, EN 14563)	60 minutes / 0.25%
Yeasts	EN 1275, EN 1650	5 minutes
	EN 13624	15 minutes
Moulds	EN 1275, T 72-300 (A. niger, A. fumigatus)	15 minutes
	EN 13624	60 minutes
Viruses	According to EN 14476 : BVDV (surrogate of HCV), PRV (surrogate of HBV), Virus Influenza A (H1N1)	5 minutes
	HIV-1	

All the product's antimicrobial activity is included in the scientific dossier, available on request.



Pavé du Moulin 59260 Lille-Hellemmes - France Tel. +33 3 20 67 67 67 - Fax : +33 3 20 67 67 68 www.anios.com

