

## Klerwipe®

# Klerwipe 70/30 Sterile Alcohol Impregnated Wipes

Range of sterile 100% knitted polyester and low particulate wipes impregnated with spore free alcohol.

### Advantages

#### ■ Rapid and effective application

The highly absorbent wipes allow a rapid clean down of surfaces, equipment and production components. The wipes remain wet throughout use ensuring a consistent application of alcohol to the whole surface area. The use of wiping can improve the effectiveness of the cleaning and disinfection procedure.



#### ■ Easy to use

The Klerwipes are folded to make them easy to remove from the pouch, even when wearing gloves. Using an impregnated wipe eliminates over-spraying of alcohol, reducing the risk of inhalation, especially in tight working spaces.

#### ■ Ideal for sterile environments

Each double bagged re-sealable pouch contains only 15 low particulate wipes, which is sufficient for most uses without compromising sterility.

#### ■ Suitable for direct product contact use

The alcohol will evaporate completely and the wipes do not contain any chemical binders, so no undesirable residues remain. This makes the Klerwipe 70/30 alcohol wipes particularly suitable for use on surfaces directly in contact with the product, where residues are unacceptable.

#### ■ Formats for all needs

There is selection of sterile alcohol impregnated wipes to suit different purposes. The pouch wipes are available in either standard 200 x 200 mm size or large 300 x 300 mm sizes. There is also a choice between either 100% polyester or low particulate wipes. The wipes are impregnated with either 70% v/v IPA or denatured ethanol blended with deionised water. In addition there are individual sachet wipes and tub wipes in various sizes.

## Description

All the different wipes are impregnated with 0.2 micron filtered alcohol, packed and double bagged in a Grade C (ISO Class 7) cleanroom, irradiated using a validated process at no less than 25kGy.

The wipes are impregnated with a choice of Klercide 70/30 IPA, a blend of 70% v/v Isopropyl Alcohol with deionised water, or Klercide 70/30 Denatured Ethanol, 70% v/v Denatured Ethanol blended with deionised water.

The 100% polyester pouch wipes are 115 gsm circular knit polyester. The low particulate wipes are made from a 68 gsm binder free, hydroentangled cellulose/polyester blend. Tub wipes are a low particulate non woven wipe.

## Pouch wipes

Each pouch contains 15 wipes so there is no wastage at the end of a manufacturing session. The polythene packaging features Klertear which makes it easy to open even when wearing gloves. The pouch is easily resealable to prevent the wipes drying out. The pouch itself does not generate any particles. 100% circular knit IPA wipes are available and also low particulate wipes, with either 70% v/v IPA or Denatured Ethanol.

The Klerwipe 70/30 alcohol pouches are also available in vaporised hydrogen peroxide impermeable packaging for use in gassing isolators.

## Sachet wipes

Individually packed low particulate wipes, folded in an easy to open foil pack. 50 foil packs are then double bagged before gamma irradiation. The foil pack eliminates lint when opening.

## Tub wipes

Two types of tub wipe are available, a high quality 42 gsm polyester/cellulose blend wipe or a 22 gsm non-woven wipe are available in a resealable tub with a low shedding aperture. Each tub is supplied double bagged and gamma irradiated. The tub contains 100 or 200 wipes respectively.

## Validation

Manufactured by a ISO 9001:2008, ISO 13485:2003 and European Medical Device Directive 93/42/EEC accredited company.

Certificates of analysis and irradiation available with every batch.

2 year shelf life.

## Product Codes

CODE	OLD CODE	DESCRIPTION	SIZE	UNIT OF SALE
3036090	WPL6500	Klerwipe 70/30 WFI Sterile 70% IPA 100% Polyester Wipe 115gsm	200x200mm	20 packs x 15 wipes
3056330		Klerwipe 70/30 WFI Sterile 70% IPA 100% Polyester Wipe 138gsm	200x200mm	20 packs x 10 wipes
3058700		Klerwipe 70/30 WFI Sterile 70% IPA 100% Polyester Wipe 138gsm	300x300mm	15 packs x 10 wipes

CODE	OLD CODE	DESCRIPTION	SIZE	UNIT OF SALE
3036040	WPK923	Klerwipe 70/30 Sterile 70% IPA Low Particulate Wipe	190x145mm	10 units x 50 sachets
3036050	WPK925	Klerwipe 70/30 Sterile 70% Denatured Ethanol Low Particulate Wipe	190x145mm	10 units x 50 sachets
3049020	WPK927	Klerwipe 70/30 Sterile 70% Denatured Ethanol Low Particulate Wipe	76x60mm	10 units x 50 sachets
3035540	WPK2500	Klerwipe 70/30 Sterile 70% IPA Low Particulate Wipe (Also available with VHP impermeable pack)	200x200mm	20 packs x 15 wipes
3035960	WPK3500	Klerwipe 70/30 Sterile 70% IPA Low Particulate Wipe	300x300mm	20 packs x 15 wipes
3035880	WPK2505	Klerwipe 70/30 Sterile 70% Denatured Ethanol Low Particulate Wipe	200x200mm	20 packs x 15 wipes
3035970	WPK3505	Klerwipe 70/30 Sterile 70% Denatured Ethanol Low Particulate Wipe	300x300mm	20 packs x 15 wipes

CODE	OLD CODE	DESCRIPTION	SIZE	UNIT OF SALE
3035980	WPK3900	Klerwipe 70/30 Sterile 70% IPA Tub Wipe 22gsm	200x230mm	15 tubs x 200 wipes
3035990	WPK3905	Klerwipe 70/30 Sterile 70% Denatured Ethanol Tub Wipe 22gsm	200x230mm	15 tubs x 200 wipes
3036000	WPK4100	Klerwipe 70/30 Sterile 70% IPA Low Particulate Tub Wipe 42gsm	200x230mm	15 tubs x 100 wipes
3036010	WPK4105	Klerwipe 70/30 Sterile 70% Denatured Ethanol Low Particulate Tub Wipe 42gsm	200x230mm	15 tubs x 100 wipes

*Use biocides safely. Always read the label and product information before use.*

*For more information please contact [info@ecolab.com](mailto:info@ecolab.com) or visit [www.ecolab.com](http://www.ecolab.com)*

