

TECHNICAL SPECIFICATION

| Klercide 70 30 Pharma Ethanol | | | |
|-------------------------------|--|--------------------|------------------|
| Product Description | Klercide 70 30 Pharma Ethanol blended with deionized water consists of a v/v formulation of 70% ethanol denatured with isopropyl alcohol, blended with 30% deionized water | | |
| Color | Colorless | | |
| Odor | Slight like alcohol | | |
| Clarity | Clear | | |
| Specific Gravity @ 20°C | 0.850 - 0.890 | | |
| Refractive Index | 19 - 21 | | |
| Production | Klercide 70 30 Pharma Ethanol blended with deionized water is filtered to 0.2 micron and bottled in a Grade C (ISO Class 7) cleanroom | | |
| Packaging | All products are cleanroom packed and double bagged. It is available in a 1 litre trigger spray fitted with the SteriShield Delivery System and in a 5 litre capped container. | | |
| Expiry | 2 years from date of manufacture | | |
| Product Codes | Code | Languages on Label | Unit of Sale |
| | 3078660 | GB; DE; FR; SLO | 6 bottles x 1.0L |
| | 3078600 | GB; DE; FR; SLO | 4 bottles x 5.0L |

USE BIOCIDES SAFELY. ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE.

PLEASE SPEAK TO YOUR ECOLAB ACCOUNT MANAGER FOR FURTHER INFORMATION OR CALL 02920 854 390 (UK) +44 2920 854 395 (EXPORT)

WORLDWIDE HEADQUARTERS

1 Ecolab Place St. Paul, MN 55102 USA www.ecolab.com/lifesciences

EUROPE HEADQUARTERS

Richtistr. 7 8304 Wallisellen Switzerland www.ecolablifesciences.com



9381.21