

The Low Residue, Broad Spectrum Rotational Disinfectant



Klercide™ Low Residue Quat Concentrate

A blend of quaternary ammonium compound and deionized water

Once appropriately diluted, Klercide Low Residue Quat Concentrate is highly efficacious at very low concentrations leaving minimal residues which are easy to remove.

Klercide Low Residue Quat Concentrate is non-corrosive and has broad spectrum activity that controls most bacterial contamination likely to be found in a cleanroom environment.

Klercide Low Residue Quat Concentrate:

- ▲ Is easy and economical to use
- ▲ Makes 20 litres of diluted disinfectant from a five litre container
- ▲ Is 0.2 micron filtered and sterile filled under unidirectional Grade A (ISO Class 5) air-flow
- ▲ Double bagged in a Grade B (ISO Class 5) cleanroom
- ▲ Is compatible with other Ecolab biocides for rotational use





Klercide™ Low Residue Quat Concentrate

Low residue, non-corrosive, broad spectrum rotational disinfectant

Why use Klercide Low Residue Quat Concentrate?

HIGH QUALITY

Klercide Low Residue Quat Concentrate is a blend of a quaternary ammonium compound and deionized water. The product is efficacious at a very low concentration ensuring the product leaves minimal residues which are easy to remove.

COST EFFECTIVE

Klercide Low Residue Quat Concentrate comes in a 5 litre container which can be diluted with the appropriate quality water to make 20 litres of diluted disinfectant. Due to being low residue, but high performance, this product can reduce the time, money and effort spent on common cleanroom issues associated with residues and their removal.

COMPATIBLE

Klercide Low Residue Quat Concentrate can be used on all types of hard surfaces in cleanrooms and controlled environments, isolators, equipment exteriors and production vessels. It is compatible with other Ecolab biocides - including Klercide Sporicidal Low Residue Peroxide and is therefore suitable for rotational use.

VALIDATION

Manufactured by a ISO 9001:2008, ISO 14001:2015 accredited company.

EN test data available on request.

A certificate of Analysis and Sterility is available for every batch. Sterility tested to an approved Pharmacopoeial method.

Product has a one year shelf life. Full supporting documentation is available for this product.



Klercide Low Residue Quat Concentrate 5L Capped

MICROBIOLOGICAL TEST RESULTS

EFFICACY UNDILUTED	CONTACT TIME
EN 1276 Bactericidal	5 min
EN 1650 Yeasticidal	15 min
EN 13697 Bactericidal and Yeasticidal	5 min

WORLDWIDE HEADQUARTERS

1 Ecolab Place
St. Paul, MN 55102
USA

www.ecolab.com

EUROPE HEADQUARTERS

Richtistr. 7
8304 Wallisellen
Switzerland

www.ecolablifesciences.com

PROVEN.
EFFECTIVE.
ECONOMICAL.

APPLICATION AREA

For use on floors, walls and hard surfaces.

INSTRUCTIONS FOR USE

Dilute before use into a clean plastic container. 5 litres should be diluted as a 1:3 dilution with water of a suitable quality to make 20 litres of diluted disinfectant. Mix well before use. For optimum results ensure complete wetting of the surfaces.

INGREDIENTS

100 ml liquid contains:

N, NDidecyl-
N, N-dimethylammoniumchloride

PACKAGING INFORMATION

UNIT SIZE	UNIT OF SALE
5.0 litre bottle	4 bottles x 5.0L

FOR FULL PRODUCT ORDER CODES, PLEASE SEE OUR PRODUCT CATALOGUE OR VISIT OUR WEBSITE

USE BIOCIDES SAFELY. ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE.

