

# The Quick and Easy Alternative to Mop and Bucket Systems



## Klerwipe™ Impregnated Sterile Mop Wipes

A unique range of impregnated mop wipes

**A unique range of sterile, impregnated mop wipes designed for use with the Klercide autoclavable mop wipe frames for optimum surface coverage.**

Impregnated with Ecolab's internationally known range of alcohols, biocides and neutral detergent, Klerwipe Sterile Impregnated Mop Wipes are a quick and convenient method of disinfecting large surfaces.

Klerwipe Sterile Impregnated Mop Wipes provide:

- ▲ A choice of wipes impregnated with, alcohols, biocides or a detergent
- ▲ Simple operation across even hard to reach areas
- ▲ Easy-open, bagged, peel-pouch packs supplied sterile
- ▲ No space consuming buckets or trolley system required
- ▲ Eliminates need to dispose of rinse and waste disinfectant solutions





# Klerwipe™ Impregnated Sterile Mop Wipes

Eliminate complex mop and bucket systems

## Why use Klerwipe Impregnated Sterile Mop Wipes?

### EFFECTIVE

All wipes are impregnated with 0.2 micron filtered fluids from the Klercide range of biocides, alcohols and detergents. Compared to bucket systems, they are:

- ▲ Faster to use
- ▲ Easier to manage
- ▲ Able to eliminate the potential risk of a dilution error when mixing a concentrate
- ▲ Provide optimal effectiveness against a wide range of commonly encountered organisms in a cleanroom

### SIMPLE

These wipes are designed to fit easily onto the Klercide Sterile Mop Wipe Frames - there are 20 wipes per pack, reducing wastage at the end of a session.

### ENVIRONMENTALLY FRIENDLY

The use of wipes rather than bucket and reservoir systems:

- ▲ Reduces wasted product
- ▲ Eliminates unnecessary mixing stages
- ▲ Managed disposal reduces the risks of biocides entering the drainage system

### VALIDATION

Manufactured by a ISO 9001:2008, ISO 14001:2015 accredited company.

Certificates of Analysis and Irradiation are available for every batch, except Klerwipe Sporidical Chlorine Quat and Sporidical Low Residue Peroxide, which are supplied with a Certificate of Sterility.

Klerwipe 70 | 30 IPA, Denatured Ethanol, Quat/Biguanide and Amine have a two year shelf life and all others have a one year shelf life.



Klerwipe Impregnated Sterile Mop Wipes

### WORLDWIDE HEADQUARTERS

1 Ecolab Place  
St. Paul, MN 55102  
USA

www.ecolab.com

### EUROPE HEADQUARTERS

Richtistr. 7  
8304 Wallisellen  
Switzerland

www.ecolablifesciences.com

EFFECTIVE.  
SIMPLE.  
CONVENIENT.

### APPLICATION AREA

Most large cleanroom surfaces. Proven effectiveness across a range of applications.

### INSTRUCTIONS FOR USE

For use with the Ecolab Klercide mop wipe frames. Pull back the upper layer of the peel-pouch pack (from left to right) and place wipe frame head on the mop wipes, wrap the wipe around the head and secure either by pushing firmly into plugs, or with stainless steel clips.

### INGREDIENTS

See individual products for details.

### PACKAGING INFORMATION

UNIT SIZE	UNIT OF SALE
420 x 250mm Wipe	10 packs x 20 wipes

FOR FULL PRODUCT ORDER CODES, PLEASE SEE OUR PRODUCT CATALOGUE OR VISIT OUR WEBSITE

USE BIOCIDES SAFELY. ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE.

