



SPECIFICALLY DESIGNED FOR FILTER MEMBRANE TESTING



Klercide™ 60|40 IPA blended with WFI

Sterile spore free 60|40 v/v IPA blended with highest grade (EP) Water for Injection

Klercide 60|40 IPA is a 60|40 v/v blend of Isopropyl alcohol and WFI is ideal for testing filter membranes and storing filters.

Filled and double bagged in a Grade C (ISO class 7) cleanroom, Klercide 60|40 IPA blended with WFI:

- ▲ Is 0.2 micron filtered
- ▲ Is spore free
- ▲ Has endotoxin levels below 0.25 EU/ml
- ▲ Is endotoxin tested to an approved pharmacopoeial method



Klercide™ 60|40 IPA

Blended with WFI



Sterile spore free 60|40 v/v IPA blended with highest grade (EP) Water for Injection

Why use the Klercide 60|40 IPA?

HIGH QUALITY

Klercide 60|40 IPA combines 60% Isopropyl alcohol blended with WFI.

The alcohol blend is 0.2 micron filtered and filled under Grade A laminar flow and double bagged in a Grade C cleanroom before irradiation using a validated process at no less than 25 kGy.

Klercide 60|40 IPA is manufactured in a facility using Quality Assurance processes to the principles of GMP to ensure reproducible quality and batch traceability. Each batch is supplied with a Certificate of Quality including sterility and irradiation.

DESIGNED FOR PURPOSE

This product has been specifically designed in line with filter manufacturers' recommendations for storage and membrane integrity testing. The bubble point test allows the detection of deviations of the required pore structure of sterile filters.

These tests are typically established by the manufacturer and demand defined solvents to achieve test results that can be compared to the filter manufacturer's tolerance data.



Klercide 60|40 IPA Blended with WFI
5 L Capped

**Safer.
Effective.
Efficient.**

APPLICATION AREA

For use in laboratory tests to measure the integrity of the filters.

INSTRUCTIONS FOR USE

Refer to filter manufacturers' instructions.

INGREDIENTS

Klercide 60|40 IPA 100ml contains: 60ml Isopropyl alcohol in Water for Injection.

PACKAGING INFORMATION

3081030 Klercide 60|40 IPA Blended with WFI

Unit Size	Unit of Sale
5 L capped	4 x 5 L

FOR PRODUCT ORDER CODES AND ADDITIONAL INFORMATION, PLEASE CONTACT YOUR ACCOUNT MANAGER OR VISIT OUR WEBSITE

USE BIOCIDES SAFELY. ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE.

WORLDWIDE HEADQUARTERS

1 Ecolab Place
St. Paul, MN 55102
USA

www.ecolab.com/lifesciences

EUROPE HEADQUARTERS

Richtistr. 7
8304 Wallisellen
Switzerland

www.ecolablifesciences.com

