



Life Sciences

USD 2956b

Cleansteam™ and Easy-Tear Cleansteam Bags



Ultra-Clean, Sterilizable and Breathable Packaging

Features

- ▶ HDPE film produced from certified virgin resins without the use of additives, in an ultra clean environment
- ▶ In-house cleaning of uncoated Tyvek 1073B
- ▶ Manufactured in ISO Class 5 cleanrooms
- ▶ Dye-penetration test (methylene blue) to show packaging integrity
- ▶ Easy-Tear option for easier opening and sealing
- ▶ Standard items available from stock; customized products on request

Benefits

- ▶ Ultra pure, clean bag with low particle, bioburden and endotoxin level (certified cleanliness level)
- ▶ Outstanding resistance to microbial penetration
- ▶ High sterility integrity until package opened
- ▶ Low particle generation when opening bag (Easy-Tear version)
- ▶ Strong, with high resistance to tear and puncture
- ▶ Sterilization indicator switching color when exposed to steam
- ▶ Increased process efficiency

Filtration. Separation. Solution.SM

Cleansteam and Easy-Tear Cleansteam bags are designed to fulfill the most stringent requirements in the sterilization process of critical components, and have outstanding performance in terms of cleanliness, material and seal strength, as well as microbial barrier properties. They are compatible with a variety of sterilization methods, such as steam and ETO sterilization, and are made from breathable, medical grade Tyvek® and own in-house extruded medical-grade HDPE.

Easy-Tear Cleansteam bags enable easy opening without the need for knives or scissors. This innovative feature virtually eliminates the risk of particle generation and accidental damage to gloves that could occur when using cutting devices, making it safer and more efficient to use in an isolator or sterile environment.

Applications

Steam or ETO sterilization of:

- ▶ Stoppers, caps and closures
- ▶ Vials and containers
- ▶ Syringe components
- ▶ Stainless steel lab or process equipment
- ▶ Medical devices
- ▶ Cleanroom garments

The bags also have an excellent chemical compatibility with Hydrogen Peroxide and are compatible with other sterilization methods such as Vaporized Hydrogen Peroxide (VPH), gamma irradiation or electron-beam.

Quality

A Certificate of Conformity (CoC) is provided with each delivery.

Regulatory compliances:

- ▶ Food contact compliant¹
- ▶ USP <87>, USP <88> Class VI and USP <661> compliant
- ▶ EP 3.1.3 compliant
- ▶ RoHS and REACH compliant
- ▶ BSE/TSE safe

¹ Manufactured from materials which are listed for food contact in FDA 21 CFR part 177 and are compliant with Commission Regulation (EU) Number 10.2011 relating to plastic materials and articles intended to come into contact with food.

Ordering Information

Cleansteam Bags

Part Number	Outer Bag Size (W x L)	Inner Bag Size (W x L)	Pkg
TH1E16P0000-0001	200 mm x 300 mm (7.9 in. x 11.8 in.)	185 mm x 250 mm (7.3 in. x 9.8 in.)	1000/ box
TH1E16P0000-0002	320 mm x 450 mm (12.6 in. x 17.7 in.)	305 mm x 400 mm (12.0 in. x 15.7 in.)	1000/ box
TH1E16P0000-0003	400 mm x 500 mm (15.7 in. x 19.7 in.)	385 mm x 450 mm (15.2 in. x 17.7 in.)	700/ box
TH1E16P0000-0004	500 mm x 650 mm (19.7 in. x 25.6 in.)	485 mm x 600 mm (19.1 in. x 23.6 in.)	500/ box
TH1E16P0000-0005	1200 mm x 1200 mm (47.2 in. x 47.2 in.)	1185 mm x 1150 mm (46.7 in. x 45.3 in.)	375/ box
TH1E16P0000-0006	650 mm x 900 mm (25.6 in. x 35.4 in.)	635 mm x 850 mm (25 in. x 33.5 in.)	250/ box

Easy-Tear Cleansteam Bags

Part Number	Outer Bag Size (W x L)	Inner Bag Size (W x L)	Pkg
TH1ET6P0000-0001	200 mm x 340 mm (7.9 in. x 13.4 in.)	185 mm x 290 mm (7.3 in. x 11.4 in.)	1000/ box
TH1ET6P0000-0002	320 mm x 450 mm (12.6 in. x 17.7 in.)	305 mm x 400 mm (12.0 in. x 15.7 in.)	1000/ box
TH1ET6P0000-0003	400 mm x 500 mm (15.7 in. x 19.7 in.)	385 mm x 450 mm (15.2 in. x 17.7 in.)	700/ box
TH1ET6P0000-0004	500 mm x 650 mm (19.7 in. x 25.6 in.)	485 mm x 600 mm (19.1 in. x 23.6 in.)	500/ box
TH1ET6P0000-0005	1200 mm x 1200 mm (47.2 in. x 47.2 in.)	1185 mm x 1150 mm (46.7 in. x 45.3 in.)	375/ box
TH1ET6P0000-0006	650 mm x 900 mm (25.6 in. x 35.4 in.)	635 mm x 850 mm (25 in. x 33.5 in.)	250/ box

Customized sizes available upon request

Width equals opening of bag

All bags are clean double packed per 100 pcs (refs -0001 to -0004), per 50 pcs (refs -0006) or per 25 pcs (refs -0005)



Life Sciences

Corporate Headquarters

Port Washington, NY, USA
+1.800.717.7255 toll free (USA)
+1.516.484.5400 phone
biopharm@pall.com e-mail

European Headquarters

Fribourg, Switzerland
+41 (0)26 350 53 00 phone
LifeSciences.EU@pall.com e-mail

Asia-Pacific Headquarters

Singapore
+65 6389 6500 phone
sgcustomerservice@pall.com e-mail

Visit us on the Web at www.pall.com/biopharm

E-mail us at allegro@pall.com

International Offices

Pall Corporation has offices and plants throughout the world in locations such as: Argentina, Australia, Austria, Belgium, Brazil, Canada, China, France, Germany, India, Indonesia, Ireland, Italy, Japan, Korea, Malaysia, Mexico, the Netherlands, New Zealand, Norway, Poland, Puerto Rico, Russia, Singapore, South Africa, Spain, Sweden, Switzerland, Taiwan, Thailand, the United Kingdom, the United States, and Venezuela. Distributors in all major industrial areas of the world. To locate the Pall office or distributor nearest you, visit www.pall.com/contact.

The information provided in this literature was reviewed for accuracy at the time of publication. Product data may be subject to change without notice. For current information consult your local Pall distributor or contact Pall directly.

© 2015, Pall Corporation. Pall, and Cleansteam are trademarks of Pall Corporation. ® indicates a trademark registered in the USA and TM is a common law trademark. **Filtration. Separation. Solution.** is a service mark of Pall Corporation. ♦Tyvek is a trademark of DuPont.