

BIOCLEAN NANO 5™ Non-sterile Nitrile Gloves NAN5

About

With textured fingers to provide enhanced grip, BioClean Nano 5™ Nitrile Gloves are latex-free, resistant to a range of chemicals and are ideal for use within electrically sensitive applications.

Specifications

COMPATIBILITY:

ISO Class 5

LENGTH:

300mm (12")

MATERIAL:

Nitrile

PROTEIN LEVEL:

No natural latex proteins

SURFACE:

Finger-textured

SHAPE:

Ambidextrous

COLOUR:

White


Features

- Resistant to a range of chemicals
- Powder-free & latex-free
- Antistatic
- Non-particulating EasyTear™ packaging

CHEMICAL RESISTANT

LATEX-FREE

ANTISTATIC



Ms.G.Love™ represents BioClean's range of cleanroom gloves

Quality Standards

- Manufactured in a facility holding ISO 9001: 2008 and 89/686/EEC Article 11b certifications
- Conforms to Category 3 - Complex Design Personal Protective Equipment PPE Directive 89/686/EEC
- Processed in an NEBB certified ISO Class 5 environment
- Complies with the requirements of Glove Standard EN420:2003 and EN374-1, 2 & 3:2003
- Physical properties comply with European Medical Glove Standard EN455-2:2009

Physical Properties

| CHARACTERISTICS | VALUE | TEST METHOD |
|--|---|---------------|
| Freedom From Holes | 1.5 AQL Performance Level 2 | EN 374-2:2003 |
| Glove Length +/-10mm (mm/inches) | 300/12 | EN 420:2003 |
| Single Wall Thickness Minimum (mils/mm) | Cuff: 3.15/0.08 Palm: 3.94/0.10 Finger Tip: 6.30/0.16 | EN 455-2:2009 |
| Force At Break During Shelf Life | ≥ 9 N | EN 455-2:2009 |

Shelf Life & Storage

Five (5) years from date of manufacture.
Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.

| PALM WIDTHS (mm) ± 5mm | | | | | | |
|------------------------|----|----|----|--------------|-----|-----|
| SIZE | XS | S | M | L | XL | XXL |
| mm | 79 | 82 | 94 | 103 | 115 | 119 |
| TEST METHOD | | | | EN455-2:2009 | | |

Glove Cleanliness Characteristics

| PARTICLES | | | | TYPICAL PARTICLE COUNT | | | | TEST METHOD |
|------------------------------------|----------|----------|----------|------------------------|---------|-----------|----------|-----------------|
| ≥ 0.5µm(counts/cm ²) | | | | <1700 | | | | IEST-RP-CC005.3 |
| TYPICAL EXTRACTABLE IONS | | | | | | | | TEST METHOD |
| ANIONS | Fluoride | Chloride | Nitrite | Bromide | Nitrate | Phosphate | Sulphate | IEST-RP-CC005.3 |
| (µg/cm ²) | ND | 0.16 | ND | ND | 0.496 | ND | 0.028 | |
| CATIONS | Lithium | Sodium | Ammonium | Potassium | Calcium | Magnesium | Zinc | |
| (µg/cm ²) | ND | 0.041 | 0.016 | 0.055 | 0.222 | ND | ND | |
| ND = Not Detected, NT = Not Tested | | | | | | | | |

TO ORDER

| Re-order code | Size |
|---------------|------|
| NAN5-XS | XS |
| NAN5-S | S |
| NAN5-M | M |
| NAN5-L | L |
| NAN5-XL | XL |
| NAN5-XXL | XXL |

Packing

100 pieces per sealed inner PE bag; one inner PE bag per sealed outer PE bag;
10 outer bags per lined carton (1000 pieces)

CONTACT

Europe (Headquarters)

Nitritex Ltd.
Minton Enterprise Park,
Oaks Drive,
Newmarket, Suffolk,
CB8 7YY,
United Kingdom
T: +44 (0) 1638 663338
F: +44 (0) 1638 668890
E: europa@bioclean.com

Other Regions

Asia: malaysia@bioclean.com
Americas: americas@bioclean.com



BioClean™ is a brand of Nitritex Ltd.

For more information, please visit www.bioclean.com

NAN5/1711/PDS9