

TECHNICAL SPECIFICATION

Klercide™ 70|30 IPA blended with WFI



Klercide 70|30 IPA blended with WFI consists of a v/v formulation of 70% isopropyl alcohol combined with 30% Water For Injection.

LIQUID SPECIFICATION	
Color	Colorless
Odor	Alcohol-like
Clarity	Clear
Endotoxin	<0.25 Endotoxin Units/ml
Specific gravity @20°C	0.872 - 0.883
IPA Content via specific gravity	68.0 - 72.2% v/v
Refractive Index @ 20°C	24 - 26

PACKAGING INFORMATION	
Production	Klercide 70 30 IPA blended with WFI is 0.2 micron filtered then filled and triple bagged in a Grade C cleanroom
Packaging	All products are cleanroom packed and triple bagged. All sprays are manufactured utilizing the SteriShield Delivery System
Sterility	Klercide 70 30 IPA blended with WFI is sterilized by gamma irradiation at not less than 25 kGy
Product Expiry	2 years from date of manufacture

PRODUCT ORDERING INFORMATION									
FORMAT	UNIT OF SALE	LANGUAGES						ORDER CODE	
500ml SDS Trigger Spray	12 x 500 ml	GB	DE	FR	IT	ES	NL	3078560	
500ml SDS Trigger Spray	12 x 500 ml	GB	SE	NO	FI	PL	DK	3078490	
1 Litre SDS Trigger Spray	6 x 1 L	GB	DE	FR	IT	ES	NL	3078540	
1 Litre SDS Trigger Spray	6 x 1 L	GB	SE	NO	FI	PL	DK	3078530	
5 Litre Capped	4 x 5 L	GB	DE	FR	IT	ES	NL	3078510	
5 Litre Capped	4 x 5 L	GB	SE	NO	FI	PL	DK	3078500	

USE BIOCIDES SAFELY. ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE.

PLEASE SPEAK TO YOUR ECOLAB ACCOUNT MANAGER FOR FURTHER INFORMATION OR CALL
02920 854 390 (UK) +44 2920 854 395 (EXPORT)

WORLDWIDE HEADQUARTERS 1 Ecolab Place, St. Paul, MN 55102, USA www.ecolab.com/lifesciences
EUROPE HEADQUARTERS Richtistr. 7, 8304 Wallisellen, Switzerland www.ecolablifesciences.com

